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Instruction Manual
Quantitative Determination of Pancreatic Elastase 1
Stool Test

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Table of Contents

	At A Glance	6
1	Introduction	21
1.1	Pathobiochemistry	21
1.2	Advantages	21
1.3	Sensitivity and specificity	22
1.4	Basic principle of the assay	23
1.5	Detection limit	23
1.6	Precision	23
1.7	Interferences	23
2	Reagents	24
3	Sample material and sample stability	24
4	Storage and stability of the test kit	25
5	Additional utensils required	25
6	Precautions	25
7	Recommendations for optimal test performance	25
7.1	Contact persons	26
8	Test procedure	27
8.1	Preparation	27
8.1.1	Preparation of sample-/washing buffer	27
8.1.2	Preparation of ELISA plate	27
8.1.3	Preparation of stool specimen	27
8.1.3.1	Using the E1 Quick-Prep™ sample preparation system	27
	Diagram of stool extraction with the E1 Quick-Prep™	28
8.1.3.2	Using the weighing method	29
8.2	Assay procedure	30
8.2.1	Incubation of samples and standards	30
8.2.2	Incubation with complex of anti-E1-bio and POD-Streptavidin	31
8.2.3	Color reaction	32
8.2.4	Stopping the color reaction	32
8.2.5	Measurement	32
8.3	Quantification of results	32
8.3.1	Manual evaluation	32
8.3.2	Evaluation by ELISA - software	33
8.3.3	Reference concentrations for pancreatic elastase 1	33
9	References	34
10	Short Protocol	Back Page

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ScheBo® • Pancreatic Elastase 1 Stool Test

Main Indications:

- Diagnosis/exclusion of exocrine **pancreatic insufficiency** caused by e.g. **Chronic Pancreatitis, Cystic Fibrosis, Diabetes Mellitus, Cholelithiasis (Gallstones), "Failure to Thrive", Pancreatic Cancer, Papillary Stenosis**
- **Follow-up monitoring** of patients with mild or moderate pancreatic insufficiency
- Diagnosis/exclusion of pancreatic involvement in association with **gastrointestinal symptoms, abdominal pain, osteoporosis**

Main Advantages:

- Digestive enzyme substitution therapy has no influence on the result of the test
- High stability allows time for convenient mailing of samples
- High specificity and sensitivity: **Specificity: 93%, Sensitivity: 93%**
- Pancreatic elastase 1 is stable during intestinal transit and absolutely pancreas-specific.
- Intra-individual variation of pancreatic elastase 1 concentration is low.

Method of Detection:

- Sandwich ELISA with two monoclonal antibodies highly specific for human pancreatic elastase 1
- The ELISA kit is based on a microtiter plate (96 well format) with **12 breakable single strips x 8 wells** suitable for up to **42 samples in duplicate**

Reference Concentration:

For adults and children after the first month of life

- Values **> 200 µg elastase/g stool** indicate normal exocrine pancreatic function
- Values **< 200 µg elastase/g stool** indicate exocrine pancreatic insufficiency

Sample Material:

- Single random stool sample (about **100 mg**) or duodenal juice
- Samples are **stable for convenient mailing** and may be stored in the laboratory up to **3 days at 4° - 8° C** or for up to **1 year at -20°C**.
- Undiluted stool extracts are stable for 1 day at 4° - 8° C or for up to 1 year at -20°C.

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1 Introduction

ScheBo® • Pancreatic Elastase 1 Stool Test for the quantitative determination of fecal levels of pancreatic elastase 1 by **healthcare professionals**.

1.1 Pathobiochemistry

Human **pancreatic elastase 1 (E1)** remains undegraded during intestinal transit. Therefore its concentration in **feces** reflects **exocrine pancreatic function**. During an inflammation of the pancreas, E1 is released into the bloodstream. Thus the quantification of pancreatic elastase 1 in **serum** allows diagnosis or exclusion of **acute pancreatitis**.

1.2 Advantages

The reliable and cost-effective ScheBo® • Pancreatic Elastase 1 Stool Test is the **new gold standard** for non-invasive pancreatic function testing. In contrast to other laboratory parameters for the diagnosis of pancreatic disease (amylase and lipase activity in serum for the diagnosis of acute pancreatitis and fecal chymotrypsin activity for the diagnosis of exocrine pancreatic insufficiency), the determination of pancreatic elastase 1 has the following advantages:

- E1 is absolutely pancreas-specific.
- Since E1 is stable during intestinal transit, the fecal elastase 1 concentration reflects the secretory capacity of the pancreas (diagnosis or exclusion of pancreatic exocrine insufficiency).
- E1 determination correlates with the gold standard invasive secretin-pancreozymin test and the secretin-caerulein test.
- Intra-individual variation of fecal E1 concentration is low.
- Enzyme substitution therapy has no influence on the determination of E1. The monoclonal antibodies used in the test do not cross-react with elastases of animal origin, which are contained in enzyme substitution preparations.
- Like other pancreatic enzymes, E1 is released into the bloodstream during pancreatic inflammation. Due to its longer half-life, compared to amylase and lipase, its concentration remains elevated longer, and an acute

pancreatitis is detectable even three or four days after onset of the disease.

Two ELISA test kits (based on monoclonal antibodies) are available for the determination of **pancreatic elastase 1**. The **serum test** quantifies E1 in serum, allowing the diagnosis or exclusion of an **acute pancreatitis** or an **inflammatory episode of chronic pancreatitis** or **ERCP- or gallstone-induced pancreatitis**. The stool test quantifies **E1 in feces** (or duodenal juice), allowing the diagnosis or exclusion of **pancreatic exocrine insufficiency**, which can be caused by **chronic pancreatitis, cystic fibrosis, diabetes mellitus, cholelithiasis (gallstones), “failure to thrive”, pancreatic cancer or papillary stenosis**, for example.

E1 serum test: cat. no.: 06

E1 stool test: cat. no.: 07

1.3 Sensitivity and specificity

The diagnostic efficiency of pancreatic elastase 1 determination in stool has been evaluated in many clinical studies. Löser *et al.* (1996) reported **93% sensitivity and 93% specificity** for the diagnosis of exocrine pancreatic insufficiency.

Stein *et al.* (1996), Löser *et al.* (1996) and Katschinski *et al.* (1997) have demonstrated a good correlation between E1 and invasive intubation tests (i.e. the secretin-pancreozymin test and the secretin-caerulein test).

In comparative studies, the ScheBo® • Pancreatic Elastase 1 Stool Test has been shown to be more sensitive than the fecal chymotrypsin test, especially in mild to moderate disease (Walkowiak *et al.* (2002), Löser *et al.* (1995 & 1996), Katschinski *et al.* (1997), Glasbrenner *et al.* (1996), Gullo *et al.* (1999)). Dominguez-Muñoz *et al.* (1995) and Glasbrenner *et al.* (1996) showed that E1 determination is more specific than the pancreolauryl test (PLT) at a comparable sensitivity. E1 has also demonstrated a better correlation with clinical response to pancreatic enzyme replacement therapy than PLT (Elphick & Kapur (2005)).

In addition, studies with E1 have shown an excellent sensitivity and specificity for the diagnosis of cystic fibrosis with pancreatic involvement (Beharry *et al.* (2002), Cade *et al.* (2000), Walkowiak *et al.* (1999 & 2002), Phillips *et al.* (1999),

Terbrack *et al.* (1996), Soldan *et al.* (1996), Gullo *et al.* (1997), Wallis *et al.* (1997)).

1.4 Basic principle of the assay

The ELISA plate is coated with a monoclonal antibody which only recognizes human pancreatic elastase 1 (E1). E1 from samples and standards binds to the antibody and is immobilized on the plate. A complex of monoclonal anti-Elastase 1-Biotin and Peroxidase (POD)-Streptavidin binds to E1 during the next incubation. The peroxidase oxidizes ABTS (2, 2'-Azino-bis-(3-ethylbenzothiazolin-6-sulfonic acid) diammonium salt), which turns dark green. Finally, the concentration of oxidized ABTS is determined photometrically.

1.5 Detection limit

Pancreatic elastase 1 is determined within the range of 15 to 500 μg E1/g stool. Concentrations below the lowest standard should be stated as $< 15 \mu\text{g}$ E1/g stool. Values above the highest standard should be indicated as $> 500 \mu\text{g}$ E1/g stool.

1.6 Precision

The intra-assay variance was evaluated by 20-fold determination of six stool samples (32-405 μg E1/g stool). The average coefficient of variance (CV) was 5.8% (range: 4.0-8.5 %).

The inter-assay variance was calculated with seven stools, which were tested on ten different days. The mean CV was 7.6 % (range: 4.8 - 10.6 %).

1.7 Interferences

Due to dilution the E1 concentration may be lowered in very watery stool samples. Therefore it is recommended to note the consistency of watery stools. In case of a pathological result ($< 200 \mu\text{g}$ E1/g stool) a formed stool sample should be requested.

There are no other factors known which interfere with the test. In particular, enzyme substitution therapy does not influence the determination of elastase 1.

2 Reagents

1. **12 breakable ELISA-strips** with 8 wells each coated with a monoclonal antibody to human pancreatic Elastase 1 (E1), 96 wells
2. **Sample-/ washing buffer concentrate (5x)** (black cap), 100 ml phosphate buffered saline, pH 7.2, with detergent
3. **E1 standards 1 to 4, ready-to-use** (blue cap), 700 µl each human pancreatic elastase 1 in aqueous solution with sodium azide
4. **Control 1, ready-to-use** (violet cap), 700 µl human pancreatic elastase 1 in aqueous solution with sodium azide
5. **Control 2, ready-to-use** (green cap), 700 µl human pancreatic elastase 1 in aqueous solution with sodium azide
6. **Complex of monoclonal anti-E1-biotin and POD-Streptavidin** (= anti-E1-bio-POD-Streptavidin-Comp.), ready-to-use, light sensitive (black plastic vial with black cap), 8 ml in aqueous solution with preservative
7. **Substrate solution, ready-to-use, light sensitive** (red cap), 12 ml ABTS in aqueous solution
8. **Stop solution, ready-to-use** (white cap), 12 ml alkaline aqueous solution

3 Sample material and sample stability

Sample material: A single random formed stool sample (a pea-sized amount is sufficient – circa 100 mg).

Sample stability: Samples are stable for convenient mailing and may then be stored in the laboratory for up to three days at 4 - 8 °C or up to one year at -20 °C. Undiluted stool extracts may be stored at 4 - 8 °C for one day or up to one year at -20 °C.

The ScheBo® • Pancreatic Elastase 1 stool test is also suitable to determine the concentration of E1 in duodenal juice (please contact us for further information). Performance characteristics have not been established for other types of samples.

4 Storage and stability of the test kit

All components of the test kit are stable at 4 - 8 °C until the expiry date shown on the kit labels. The kit must not be used after the expiry date. Unused ELISA-strips and/or wells must be stored in the re-sealable foil bag containing the desiccant, making sure the seal is completely closed.

5 Additional utensils required

- Polystyrene test tubes (5 ml, 10 ml and 12 ml) with caps
- 500 ml graduated cylinder
- vortex mixer
- adjustable precision pipettes: 0-50 µl, 50-200 µl, and 200-1000 µl
- 5 ml and 10 ml pipettes
- adjustable 8 channel pipette 50-250 µl
- ELISA reader capable of reading absorbance at 405 nm.

6 Precautions

For *in vitro* diagnostic use only. Extraction buffer, standards, control and the complex of anti-E1-bio and POD-Streptavidin contain a preservative. **Please observe the safety regulations. Avoid skin contact. Do not pipette by mouth. Wear disposable gloves** while performing the tests. A materials safety data sheet is available on request. Do not mix materials from different master lots.

7 Recommendations for optimal test performance

1. Follow the instruction manual precisely.
2. All components of the test kit have to be stored at 4 - 8 °C. Bring all reagents including the ELISA plate to room temperature shortly before use. After using store all reagents at 4 - 8 °C immediately.
3. Vortex all liquid reagents before use. Avoid droplets in the caps of the tubes.
4. Do not touch the bottom of the ELISA plate.
5. Always work up the ELISA plate in identical order and equal time intervals.

6. To avoid contamination use clean pipette tips and clean receptacles. Do not use the same receptacle for different reagents (e.g substrate and stop solution).
7. Pipetting
 - Pipette all reagents and samples into the lower third of a well
 - For washing hold the pipette tips at the upper rim of the wells
8. Washing
 - Before each washing step, invert the plate and **tap hard** on a clean paper towel to **remove any remaining liquid**.
 - Use only clean paper towels for draining off liquid.
 - Make sure that no liquid drains back.
 - Incubate with washing buffer at least **1-2 minutes per washing step**.
9. Measurement
 - Agitate plate well before each measurement to ensure even distribution of the dye.
 - **Remove any bubbles** with a clean needle.
 - Wait at least five minutes after adding the stop solution before measuring.

7.1 Contact persons

Dr. Karin Decker

ScheBo® • Biotech AG

Netanyastr. 3, 35394 Giessen, Germany

Tel: +49-(0)641-4996-0, Fax: +49-(0)641-4996-77

<http://www.schebo.de>, e-mail: schebo@schebo.com

- your local distributor

8 Test procedure

8.1 Preparation

8.1.1 Preparation of sample-/washing buffer

100 ml sample-/washing buffer 5x (black cap) + 400 ml bidistilled water.
The diluted sample-/washing buffer is stable for 6 months at 4 - 8 °C.

8.1.2 Preparation of ELISA plate

Bring ELISA plate to room temperature before opening. Tear open the resealable foil bag where indicated by the notch, open the re-sealable closure and remove desired number of breakable ELISA strips. **Unused ELISA strips and/or wells must be stored in the resealable foil bag containing the desiccant, making sure the seal is completely closed.**

8.1.3 Preparation of stool specimen

- We recommend the E1 Quick-Prep™ dosing device (cat. no. 07-Quick, see 8.1.3.1) is used for speed and convenience, or alternatively
- Weighing method: stool specimen can be weighed (see 8.1.3.2)

8.1.3.1 Using the E1 Quick-Prep™ sample preparation system

The tubes of the ScheBo® • E1 Quick-Prep™ sample preparation system contain ready-to-use extraction buffer. See the imprint on the packaging for the expiry date.

Preparation of stool specimens

The ScheBo® • E1 Quick-Prep™ (cat. no. 07-Quick) should be used according to figure 1 (page 28).

Dilution of stool extracts (1:70)

Preparation of 1:70 dilution:

10 µl extracted stool sample+ 700 µl sample-/washing buffer

Stool extraction with E1 Quick-Prep™


Figure 1

8.1.3.2 Performance with the weighing method

Extraction buffer concentrate (5x) for stool specimen (square flask, green cap), 100 ml (**cat. no.: 02**) phosphate buffered saline, pH 7.2, with detergent and sodium azide

Preparation of extraction buffer

Dilute the extraction buffer (5x) (cat. no.: 02): 100 ml extraction buffer 5x + 400 ml bidistilled water. The diluted extraction buffer is stable for 6 months at 4 - 8 °C.

Weighing the stool specimen

Tare tube and inoculating loop on a sensitive digital laboratory balance. Take a small sample (approximately 100 mg) from the stool specimen with the inoculating loop and replace the loop into the tube to weigh the sample. A wooden applicator or tooth pick may be used instead of the inoculating loop. Add extraction buffer to the stool sample according to the mass of the sample, (e.g. 100 mg stool + 10 ml extraction buffer or 75 mg stool + 7.5 ml extraction buffer). The final concentration must be 10 mg stool/ml extraction buffer.

Homogenisation and extraction of stool samples

Vortex the stool suspension thoroughly at room temperature. Stool suspensions must be homogenized to ensure a complete extraction of pancreatic elastase 1. After an extraction period of at least five minutes, vortex the suspension again. Then, after any particles have settled, dilute the stool extracts. Because pancreatic elastase 1 is very stable, the extraction period may be extended for up to 24h at 4 - 8 °C and the dilutions can be performed the next day.

Dilution of stool extracts (1:250)

Preparation of 1:250 dilution:

10 µl extracted stool sample + 2.5 ml sample-/washing buffer.

8.2 Assay procedure

8.2.1 Incubation of samples and standards

	1	2	3	4	5	6	7	8	9	10	11	12
A	Blank	Blank	S3	S3	S11	S11	S19	S19	S27	S27	S35	S35
B	STD1	STD1	S4	S4	S12	S12	S20	S20	S28	S28	S36	S36
C	STD2	STD2	S5	S5	S13	S13	S21	S21	S29	S29	S37	S37
D	STD3	STD3	S6	S6	S14	S14	S22	S22	S30	S30	S38	S38
E	STD4	STD4	S7	S7	S15	S15	S23	S23	S31	S31	S39	S39
F	CON	CON	S8	S8	S16	S16	S24	S24	S32	S32	S40	S40
G	S1	S1	S9	S9	S17	S17	S25	S25	S33	S33	S41	S41
H	S2	S2	S10	S10	S18	S18	S26	S26	S34	S34	S42	S42

Figure 2: Possible plate layout

STD: standards

CON: control

S1-S42: samples

Blank = wells A1 and A2; pipette 50 μ l of sample-/washing buffer into each well.

Standards (blue cap) are ready-to-use; pipette 50 μ l of each standard (undiluted) into columns 1 and 2 as duplicate.

Stool

Standard 1 corresponds to 15 μ g/g

Standard 2 corresponds to 50 μ g/g

Standard 3 corresponds to 150 μ g/g

Standard 4 corresponds to 500 μ g/g

Control 1, ready-to-use (violet cap); pipette 50 μ l into wells F1 and F2. Control 1 must be determined on each run.

Control 2, ready-to-use (green cap).

Control 2 may be used if required by local, state or national regulations/guidelines or accrediting organisations.

Pipette 50 µl into wells G1 and G2 (instead of a sample).

Stool

Control 1 corresponds to 200 µg/g ± 15 %

Control 2 - see label for concentration

Pipette 50 µl of extracted diluted stool specimen (see 8.1.3.1 und 8.1.3.2) from each sample into each of two adjacent wells. The control can also be used as an internal control.

Incubate for 30 minutes at room temperature.

Washing: Empty the wells and wash each well 3 times with sample-/washing buffer (8 channel pipette, 250 µl/well). Invert the plate and **tap it hard** on a clean paper towel to **completely remove any remaining liquid**.

8.2.2 Incubation with complex of anti-E1-bio and POD-Streptavidin

Add 50 µl/well ready-to-use complex of anti-E1-biotin and POD-Streptavidin = anti-E1-bio-POD-Streptavidin-Comp. (black plastic vial with black cap)

Incubate for 15 minutes in the dark at room temperature.

Washing: Empty the wells and wash each well 3 times with sample-/washing buffer (8 channel pipette, 250 µl/well). Invert the plate and **tap it hard** on a clean paper towel to **completely remove any remaining liquid**.

8.2.3 Color Reaction

Add 100 µl of ready-to-use substrate solution (red cap) to each well.

Incubate for 15 minutes in the dark at room temperature.

(You may need to shorten this time when using an ELISA plate-reader or a fully automated machine which reads absorbances up to only 2.5 or 3.)

8.2.4 Stopping the color reaction

Stop the substrate reaction by adding 100 µl of stop solution per well (ready-to-use, white cap). Mix contents well by agitating the plate.

8.2.5 Measurement

Read the optical density at **405 nm** with a microtiter plate reader **between 5 and 30 minutes** after addition of the stop solution. Mix contents well before measuring. **492 nm** can be used as a reference wavelength.

8.3 Quantification of results

8.3.1 Manual evaluation

Calculate the mean optical densities of all duplicates after mean blank value has been subtracted. Plot the concentration of standards versus their corresponding optical densities on a log-log paper = standard curve.


Figure 3: Typical example of a standard curve

Example: The average optical density after the mean blank value has been subtracted was 0.65. This optical density corresponds to a concentration of 120 µg E1/g stool (see standard curve).

8.3.2 Evaluation by ELISA - software

Define blank, standards and samples according to the plate layout (figure 2). Use the curve - fit method (linear regression) with log - log scale.

8.3.3 Reference concentrations for pancreatic elastase 1

In stool:

normal:	200 to >500 µg E1/g stool
moderate to mild exocrine pancreatic insufficiency:	100 – 200 µg E1/g stool
severe exocrine pancreatic insufficiency:	<100 µg E1/g stool

These pancreatic elastase 1 concentrations refer to formed stool samples. In case of pathological elastase 1 concentrations (< 200 µg E1/ g stool) in watery stool samples a second formed stool sample should be requested (see section 1.7).

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- Hardt, P.D., Marzeion, A.M., Schnell-Kretschmer, H., Wüsten, O., Nalop, J., Zekorn, T., Klör, H.U. (2002) Fecal elastase 1 measurement compared with endoscopic retrograde cholangiopancreatography for the diagnosis of chronic pancreatitis, *Pancreas* 25: e6-e9
- Henker, J., Thamm, K. (1999) Untersuchungen zur exokrinen Pankreasfunktion bei Kindern und Jugendlichen mit Diabetes mellitus Typ 1, Abstract/ 5. Ärzteswoche Gastroenterologie/Viszeralchirurgie in Garmisch-Partenkirchen, 21.-26.02.1999
- Icks, A, Haastert, B., Giani, G., W. Rathmann (2001) Low fecal elastase-1 in type I diabetes mellitus, *Z Gastroenterol* 39: 823-830
- Katschinski, M., Schirra, J., Bross, A., Göke, B., Arnold, R. (1997) Duodenal Secretion and Fecal Excretion of Pancreatic Elastase-1 in Healthy Humans and Patients with Chronic Pancreatitis, *Pancreas* 15,2: 191-200
- Leeds, J. S., Hopper, A. D., Hurlstone, D. P., Edwards, S. J., McAlindon, M. E., Lobo, A. J., Donnelly, M. T., Morley, S., Sanders, D.S. (2007) Is exocrine pancreatic insufficiency in adult coeliac disease a cause of persisting symptoms? *Aliment Pharmacol Ther* 25, 265–271
- Leus, J., Van Biervliet, St., Robberecht, E. (2000) Detection and Follow up of Exocrine Pancreatic Insufficiency in Cystic Fibrosis: A Review, *Eur J Pediatr* 159: 563-568
- Löser, Chr., Brauer, C., Aygen, S., Hennemann, O., Fölsch, U.R. (1998) Comparative Clinical Evaluation of the ¹³C-Mixed Triglyceride Breath Test as an Indirect Pancreatic Function Test, *Scand J Gastroenterol* 33: 327-334
- Löser C., Fölsch U. R. (1996) Diagnostik der chronischen Pankreatitis, *Dtsch med Wschr* 121: 243-247
- Löser C., Mølgaard, A., Fölsch, U. R. (1996) Faecal elastase 1: a novel, highly sensitive, and specific pancreatic function test, *Gut* 39, 580-586
- Malfertheiner, P., Büchler, M., Stanescu, A., Uhl, W., Ditschuneit, H. (1987) Serum elastase 1 in inflammatory pancreatic and gastrointestinal diseases and in renal insufficiency - A comparison with other serum pancreatic enzymes, *Int J Pancreat* 2: 159 - 170

- Nissler, K., Von Katte, I., Huebner, A., Henker, J. (2001) Pancreatic Elastase 1 in Feces of Preterm and Term Infants, *J Ped Gastroenterol Nutr* 33: 28-31
- Nunes, A. C .R., Pontes, J.M., Rosa, A., Gomes, L., Carvalheiro, M., Freitas, D. (2003) Screening for Pancreatic Exocrine Insufficiency in Patients with Diabetes Mellitus, *Am J Gastroenterol* 98 (12): 2672-2675
- Peters, K., Lerche, L., Stange, E.F. (1998) Erhöhte Serum-Elastase als Risikofaktor für eine Post-ERCP-Pankreatitis *Z Gastroenterol.* 36: 697
- Phillips, I.J., Rowe, D.J., Connett, G., Dewar, P. (1999) Faecal elastase 1: A marker of exocrine pancreatic insufficiency in cystic fibrosis, *Ann Clin Biochem* 36:739-742
- Rathmann, W., Haastert, B., Icks, A., Giani, G., Hennings, S., Mitchell, J., Curran, S., Wareham, N.J., (2001) Low faecal elastase 1 concentrations in type 2 diabetes mellitus, *Scand J Gastroenterol* 36 (10):1056-61
- Schappi, M.G., Smith, V.V., Cubitt, D., Milla, P.J., Lindley K.J. (2002) Faecal Elastase 1 Concentration is a Marker of Duodenal Enteropathy, *Arch Dis Child* 86(1):50-53
- Scheefers-Borchel, U., Scheefers, H., Arnold, R., Fischer, P., Sziegoleit, A. (1992) Pankreatische Elastase 1: Parameter für die chronische und akute Pankreatitis, *Lab. med.* 16: 427 – 432
- Salvatore, S., Finazzi, S., Barassi, A., Verzelletti, M., Tosi, A., d'Eril, G. V.M., Nespoli, L. (2003) Low Fecal Elastase: Potentially Related to Transient Small Bowel Damage Resulting from Enteric Pathogens, *J Ped Gastroenterol Nutr* 36: 392-396
- Seno, T., Harada, H., Ochi, K., Tanaka, J., Matsumoto, S., Choudhury, R., Mizushima, T., Tsuboi, K. & Ishida, M. (1995) Serum Levels of Six Pancreatic Enzymes as Related to the Degree of Renal Dysfunction, *Am J Gastroenterol* 90: 2002-2005
- Soldan, W., Henker, J., Sprössig, C. (1997) Sensitivity and Specificity of quantitative Determination of Pancreatic Elastase 1 in Feces of Children, *J Ped Gastroenterol Nutr* 24: 53-55
- Stein, J. (1996) Where to place Human Elastase-1 in the Diagnosis of Pancreatic Insufficiency, *Z Gastroenterol.* 34 (S4): 6-10

- Stein, J., Jung, M., Sziegoleit, A., Zeuzem, S., Caspary, W.F., Lembcke, B. (1996) Fecal immunoreactive elastase-1: Evaluation of a new tubeless pancreatic function test in comparison with other indirect and direct tests for exocrine pancreatic function Clin Chem 42: 222-226
- Stein, J., Spirchez, Z., Lembcke, B., Caspary, F. (1997) Untersuchungen zur Bedeutung der Elastase-Bestimmung als einen neuen nicht-invasiven Test der exokrinen Pankreasinsuffizienz, Z Gastroenterol. (Suppl. 1): 122-129
- Stevens, T., Conwell, D., Zuccaro, G., van Lente, F., Khandwala, F., Hanaway, P., Vargo, J.J., Dumot, J.A. Analysis of pancreatic elastase-1 concentrations in duodenal aspirates from healthy subjects and patients with chronic pancreatitis, Dig Dis Sci 49: 1405-1411
- Sziegoleit, A., Krause, E., Klör, H.U., Kranacher, L., Linder, D. (1989) Elastase 1 and Chymotrypsin B in pancreatic juice and feces, Clin Biochem. 22: 85 - 89
- Sziegoleit, A., Linder, D. (1991) Die menschliche Pankreaselastase 1, medwelt 42: 682 - 684
- Teichmann, J., Lange, U., Hardt, P.D., Schnell-Kretschmer, H., Stracke, H., Bretzel, R.G., Klör, H.U. (2001) Pancreatic Elastase 1 in Patients with Osteoporosis, Abstract Congress Osteologie 2001, 4-7. April 2001, Wiesbaden
- Terbrack, H.G., Gürtler, K.H., Hüls, G., Bittner-Dersch, P., Klör, H.U., Lindemann, H. (1996) Humanspezifische fäkale Pankreaselastase bei Kindern, Monatsschr. Kinderheilkd. 144: 901-905
- Terbrack, H.G., Gürtler, K.H., Klör, H.U., Lindemann, H. (1995) Human Pancreatic Elastase 1 Concentration in Faeces of Healthy Children and Children with Cystic Fibrosis, Gut, 4th UEGW Berlin 2262: A 253
- Walkowiak, J. (2004) Assessment of maldigestion in cystic fibrosis, J Pediatr 145: 285-287
- Walkowiak, J., Cichy, W.K., Herzig, K.H. (1999) Comparison of Fecal Elastase-1 Determination with the Secretin-Cholecystokinin Test in Patients with Cystic Fibrosis Scand, J Gastroenterol 34(2): 201-207

- Walkowiak, J. (2000) Faecal Elastase 1: Clinical Value in the Assessment of Exocrine Pancreatic Function in Children, *Eur J Ped* 159: 869-870
- Walkowiak, J., Herzig, K.-H. (2001) Fecal Elastase-1 is decreased in villous atrophy regardless of the underlying disease, *Eur J Clin Invest* 31:425-430
- Walkowiak, J., Herzig, K.-H., Strzykala, K., Przyslawski, J. and Krawczynski, M. (2002) Fecal Elastase-1 Is Superior to Fecal Chymotrypsin in the Assessment of Pancreatic Involvement in Cystic Fibrosis, *Pediatrics* 110: 1-4
- Walkowiak, J., Nousia-Arvanitakis, S., Henker, J., Stern, M., Sinaasappel, M., Dodge J.A. (2005) Indirect Pancreatic Function Tests in Children, *J Pediatr Gastroenterol Nutr* 40: 107–114
- Wallis, C., Leung, T., Cubitt, D., Reynolds, A. (1997) Stool elastase as a diagnostic test for pancreatic function in children with cystic fibrosis, *The Lancet* 350: 1001

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Short protocol for the experienced user

Important: The short protocol is not a substitute for the detailed protocol given in this leaflet !

ScheBo® • Pancreatic Elastase 1 Stool Test for the quantitative determination of fecal levels of pancreatic elastase 1 by **healthcare professionals**.

- Prepare the sample-/washing buffer and the extraction buffer
- Extract and homogenize stool - Dilute stool extract in sample-/washing buffer
- Pipette 50 µl blank, standards, control and samples in duplicate in the ELISA-strips - Incubate 30 minutes at room temperature - Wash
- 50 µl anti-E1-bio and POD-Streptavidin-Complex (ready-to-use) - Incubate 15 minutes at room temperature (in the dark) - Wash
- 100 µl substrate solution (ready-to-use) - Incubate 15 minutes at room temperature (in the dark)
- Add 100 µl stop solution (ready-to-use)
- Read plate at OD_{405} or $OD_{405} - OD_{492}$
- Evaluate with standard curve using a log-log scale